

# SURGICAL HEARING AIDS

BY,  
G SUDHANTHIRA M.E.,  
A/P BME

# Comparison of Surgical vs. Non-Surgical Hearing Aids

Feature	Non-Surgical Hearing Aids	Surgical Hearing Aids
Treatment Type	External devices	Implantable devices
Indications	Mild to severe hearing loss	Severe to profound hearing loss
Customization	Programmable settings	Requires surgery but permanent
Suitability	Most individuals	Limited to specific hearing conditions
Effectiveness	Amplifies sound	Restores hearing in cases where hearing aids fail
Cost	Affordable to moderate cost	Expensive (requires surgery)
Maintenance	Requires battery change and cleaning	Regular follow-ups needed
Cosmetic Appeal	Visible in some cases	Invisible (fully implantable)

# Surgical Hearing Aids

- Surgical hearing aids, also known as implantable hearing devices, are devices that require surgery to be implanted in the ear or skull for individuals with severe to profound hearing loss or certain ear conditions. Surgical hearing aids bypass the outer and middle ear, offering a more permanent solution to hearing loss.
- **Types:**
  - Cochlear Implants (CI)
  - Bone-Anchored Hearing Aids (BAHA)
  - Middle Ear Implants (MEI)
  - Auditory brainstem implants (ABI)

# Cochlear Implants



- The cochlea is a spiral-shaped, fluid-filled structure in the inner ear responsible for converting sound waves into electrical impulses that are interpreted by the brain. This process, known as auditory
- Cochlear implants (CI) are advanced medical devices designed to provide hearing to individuals with **severe to profound sensorineural hearing loss** who do not benefit from traditional hearing aids. Unlike hearing aids, which amplify sound, cochlear implants bypass the damaged parts of the inner ear and directly **stimulate the auditory nerve** using electrical signals.

# Structure of a Cochlear Implant

A cochlear implant consists of **two main components**:

1. **External Components** (Worn outside the body):
  - o Microphone – Captures sound from the environment.
  - o Speech Processor – Converts sound into digital signals.
  - o Transmitter Coil – Sends the processed sound signals to the internal implant via radio waves.
2. **Internal Components** (Surgically implanted):
  - o Receiver-Stimulator – Placed under the skin behind the ear, it receives signals from the transmitter.
  - o Electrode Array – Inserted into the cochlea (inner ear), it stimulates the auditory nerve directly, sending sound signals to the brain.

# Cochlear Implants



# Mechanism:

## Step 1: Sound Capture and Processing

- The microphone picks up external sound and sends it to the speech processor.
- The speech processor filters, analyzes, and converts the sound into digital signals using advanced processing algorithms.
- The processed signals are transmitted via radio waves to the internal receiver-stimulator through the transmitter coil.

## Step 2: Signal Transmission to the Cochlea

- The receiver-stimulator converts the digital signals into electrical impulses.
- These impulses travel through the electrode array implanted in the cochlea.
- The electrodes stimulate different regions of the auditory nerve based on the tonotopic organization of the cochlea (high frequencies at the base, low frequencies at the apex).

# **Mechanism: Cochlear Implants**

## **Step 3: Auditory Nerve Activation**

- The electrical impulses directly stimulate the spiral ganglion neurons of the auditory nerve.
- The auditory nerve fibers carry these signals to the brainstem and eventually to the auditory cortex in the brain.
- The brain interprets these signals as perceived sound.

## **Role of Frequency and Loudness Encoding**

- Frequency Encoding: Different electrodes along the electrode array stimulate different cochlear regions, mimicking the natural frequency mapping (tonotopy) of the cochlea.
- Loudness Encoding: The strength of electrical stimulation determines the perceived loudness of the sound.



# Cochlear Implants

## **Advantages:**

- Provides hearing ability to people with profound deafness.
- Helps children develop speech and language skills.
- Improves communication and social interaction.
- Enhances sound localization and awareness.
- Long-term solution compared to hearing aids.

## **Limitations & Risks:**

- Does not restore normal hearing; it provides a different way of hearing.
- Requires surgery, which carries risks like infection and bleeding.
- Sounds may initially seem mechanical or artificial.
- Requires extensive rehabilitation for the best results.
- Costly procedure, and maintenance is required.

# Cochlear Implants

Feature	Hearing Aid	Cochlear Implant
<b>Function</b>	Amplifies sound for better hearing	Bypasses damaged ear structures and stimulates the auditory nerve directly
<b>Best for</b>	People with mild to moderate hearing loss	People with severe to profound hearing loss
<b>Working Principle</b>	Amplifies external sounds and sends them to the ear	Converts sound into electrical signals and sends them to the brain
<b>External Components</b>	Microphone, amplifier, speaker	Microphone, speech processor, transmitter
<b>Internal Components</b>	None (fully external device)	Receiver and electrode array implanted in the cochlea
<b>Surgical Requirement</b>	No surgery needed	Requires surgery to implant internal components
<b>Reversibility</b>	Can be removed anytime	Permanent, but can be turned off
<b>Cost</b>	Generally lower cost	Higher cost due to surgery and device complexity
<b>Suitability for Deaf Individuals</b>	Limited benefit for profound hearing loss	More effective for profound hearing loss
<b>Sound Quality</b>	Natural sound but may be distorted for severe loss	Different from natural hearing; brain adapts over time
<b>Maintenance</b>	Regular battery or recharge	Requires both external and internal maintenance

# Bone-Anchored Hearing Aids (BAHA)

Bone-Anchored Hearing Aids (BAHA) are implantable hearing devices designed for individuals with conductive, mixed hearing loss, or single-sided deafness (SSD). Unlike traditional hearing aids that amplify sound through the ear canal, BAHA uses bone conduction to transmit sound directly to the inner ear.

## **Structure:**

A BAHA system consists of three main components:

### **1.External Sound Processor**

- Captures sound from the environment.
- Converts sound into vibrations.
- Sends these vibrations to the titanium implant.

### **2.Titanium Implant (Screw)**

- Surgically placed into the skull bone (mastoid).
- Acts as an anchor for the sound processor.
- Helps transmit vibrations to the cochlea.

### **3.Abutment or Magnetic Attachment**

- Connects the external sound processor to the titanium implant.
- Can be a direct connection (abutment) or a magnetic coupling under the skin.

# Types of BAHA Attachments


## **1. Abutment-Based BAHA (Percutaneous)**

- A small abutment sticks out through the skin.
- The sound processor directly attaches to it.
- Provides better sound transmission but requires regular cleaning.


## **2. Magnetic BAHA (Transcutaneous)**

- The implant is fully under the skin.
- The sound processor attaches via a magnet.
- More cosmetic and comfortable, but sound quality may be slightly reduced

# Transcutaneous



# Percutaneous



# Mechanism of Percutaneous BAHA

## Step 1: Sound Capture and Processing

- The microphone picks up environmental sounds.
- The sound processor amplifies and converts the sounds into mechanical vibrations.
- The vibrations are **transmitted through a titanium abutment**, which directly connects the processor to the skull bone.

## Step 2: Direct Bone Conduction

- The titanium implant transmits vibrations directly into the skull bone without any soft tissue interference.
- These vibrations reach the cochlea, bypassing the outer and middle ear.

## Step 3: Cochlear Stimulation and Hearing Perception

- The cochlea detects vibrations, and hair cells convert them into electrical signals.
- These signals are transmitted through the auditory nerve to the brain, where they are interpreted as sound.

# Mechanism of Transcutaneous BAHA

## Step 1: Sound Capture and Processing

- The microphone captures external sounds.
- The sound processor amplifies and converts sounds into mechanical vibrations.
- The vibrations are transmitted through **magnetic coupling or a passive connection**, eliminating the need for a skin-penetrating abutment.

## Step 2: Indirect Bone Conduction

- The internal implant receives the vibrations and transfers them to the mastoid bone through soft tissue.
- The vibrations propagate through the skull bone to the cochlea, bypassing the outer and middle ear.

## Step 3: Cochlear Stimulation and Hearing Perception

- The cochlea perceives the vibrations as sound, just like in natural hearing.
- Hair cells convert these mechanical vibrations into neural signals.
- The signals travel via the auditory nerve to the brain, allowing sound perception.

# Comparison: Transcutaneous vs. Percutaneous BAHA

Feature	Percutaneous BAHA	Transcutaneous BAHA
Connection Type	Direct abutment (skin-penetrating)	Magnetic coupling (skin intact)
Sound Transmission	Direct bone conduction (stronger)	Indirect bone conduction (weaker)
Surgical Complexity	Requires abutment attachment	Requires internal magnet implant
Risk of Infection	Higher (skin penetration)	Lower (no skin breach)
Cosmetic Appeal	Visible abutment	More discreet
Ideal for	Severe conductive/mixed hearing loss	Mild-to-moderate conductive hearing loss



# BAHA Systems

## **Advantages of BAHA Systems:**

- Bypasses the damaged outer and middle ear, making it ideal for conductive hearing loss.
- Provides natural sound perception through bone conduction.
- Works for single-sided deafness (SSD) by transferring sound to the healthy ear.
- Improved speech clarity in noisy environments compared to conventional hearing aids.

## **Limitations of BAHA Systems:**

- Surgical Risks: Infection, implant failure, or skin irritation.
- Sound Quality Difference: Percutaneous BAHA provides clearer sound than transcutaneous BAHA.
- High Cost: Expensive compared to traditional hearing aids.
- Requires Maintenance: External components need regular care.

## **Applications of BAHA:**

Conductive Hearing Loss: Due to chronic ear infections, atresia, or malformations.

Mixed Hearing Loss: When both conductive and sensorineural components are present.

Single-Sided Deafness (SSD): Transmits sound from the deaf side to the functional cochlea.

# Middle Ear Implants (MEI)

An MEI is a surgically implanted device that works by **vibrating the bones** of the middle ear (ossicles) or the cochlea directly to transmit sound. This is particularly useful for individuals with sensorineural hearing loss (damage to the cochlea or auditory nerve) or conductive hearing loss (problems in the outer or middle ear). The device enhances sound by amplifying vibrations and sending them to the cochlea for auditory perception.

## **Advantage of MEIs:**

- Do not obstruct the ear canal.
- No occlusion effect.
- Clearer sound quality.
- More stable and reliable

# Types of Middle Ear Implants


## **1. Fully Implantable MEIs**

- The Esteem Implant is a fully implantable device that is placed entirely inside the ear, with no external parts visible. It consists of a microphone, sound processor, receiver, and vibrating component that are all embedded in the ear.
- The microphone captures sound and sends it to the sound processor, which amplifies it and converts it into a signal that stimulates the eardrum. The eardrum vibrations are transferred to the ossicles.

## **2. Semi-Implantable MEIs:**

- The Vibrant Soundbridge is a semi-implantable device that has an external sound processor attached to a receiver/transducer implanted in the middle ear. The receiver is usually placed on the ossicles (the bones in the middle ear) or close to the cochlea.
- The device amplifies sound and sends mechanical vibrations through the bones of the middle ear, bypassing damaged parts of the ear, such as the eardrum or ossicles.

# Middle Ear Implants (MEI)



# Components of Middle Ear Implants (MEI)

MEIs are typically made up of two main parts:

## 1. External Components:

- **Microphone:** The microphone is the external component that picks up sound from the environment. It is usually attached to a small processor, typically worn behind the ear, similar to a hearing aid.
- **Sound Processor:** The sound processor is the part of the external component that converts the **captured sound into an electrical signal**. This signal is transmitted to the internal device through **a magnetic connection or transcutaneous link**.

## 2. Internal Components:

- **Receiver/Stimulator:** The receiver is implanted inside the ear, often placed beneath the skin behind the ear. It receives the electrical signals from the external processor and converts them into mechanical vibrations.
- **Vibrating Element (Vibratory Transducer):** The internal vibrating component directly stimulates the ossicles (the small bones in the middle ear) or attaches to the cochlea. This component transforms the electrical signals from the receiver into physical vibrations that are then transmitted to the cochlea.

# Auditory Brainstem Implant (ABI)

- An **Auditory Brainstem Implant (ABI)** is a specialized neuroprosthetic device designed for individuals who cannot benefit from cochlear implants due to damage or absence of the auditory nerve. This device is primarily used in patients with **Neurofibromatosis Type 2 (NF2)**, congenital inner ear malformations, or auditory nerve damage.
- **Components:**
  1. **External Components (Similar to a cochlear implant)**
    - **Microphone:** Captures sound from the environment.
    - **Speech Processor:** Converts sound into electrical signals.
    - **Transmitter Coil:** Sends processed signals wirelessly to the internal implant.
  2. **Internal Components**
    - **Receiver-Stimulator:** Implanted under the skin to receive signals from the external processor.
    - **Electrode Array:** Placed directly on the **cochlear nucleus** in the brainstem to stimulate neural activity.

# Mechanism of Action

- The **microphone** picks up sound and sends it to the **speech processor**.
- The **speech processor** converts sound into electrical signals and transmits them to the **receiver-stimulator** via the **transmitter coil**.
- The **receiver-stimulator** sends electrical impulses to the **electrode array** implanted on the **cochlear nucleus**.
- These impulses stimulate the neurons in the cochlear nucleus, allowing the brain to interpret sound.



# Auditory Brainstem Implant (ABI)

## **Advantages of Auditory Brainstem Implants**

- Bypasses the auditory nerve and directly stimulates the cochlear nucleus.
- Useful for patients who cannot benefit from cochlear implants.
- Helps improve sound awareness, speech perception, and communication skills.
- Can be implanted in both children and adults.

## **Limitations and Challenges**

- Surgical risks: Infection, cerebrospinal fluid (CSF) leakage, facial nerve injury.
- Limited speech perception: ABI provides sound awareness but not natural hearing.
- Expensive and less available than cochlear implants.
- Requires extensive rehabilitation and training.