

Devices to Improve Hearing and Assistive Listening Devices (ALDs)

BY,
G SUDHANTHIRA M.E.,
A/P BME

Definition:

- **Devices to Improve Hearing** are primarily medical or technological devices designed to enhance an individual's hearing abilities by amplifying sounds and improving clarity. These devices are primarily used by people who suffer from hearing impairments, and they focus on amplifying sounds so that individuals can hear more effectively in everyday situations.
- **Assistive Listening Devices (ALDs)**, on the other hand, are specialized systems that help improve sound quality in specific listening environments. These devices are often used in situations where background noise is present, or the distance between the listener and the sound source is too great for standard hearing aids to be effective. ALDs can be used alongside hearing aids or cochlear implants, or they can function independently to improve the quality of speech comprehension, particularly in challenging auditory environments.

Function and Purpose

- The primary function of **Devices to Improve Hearing** is to compensate for hearing loss by amplifying sound. These devices are typically customized to meet the individual's specific hearing profile, making them an ideal solution for people with mild to profound hearing impairments. Their focus is on improving overall hearing ability and sound clarity
- In contrast, the main function of **Assistive Listening Devices (ALDs)** is to improve sound clarity in specific, often noisy, environments. These devices are ideal for settings such as lecture halls, movie theaters, or places of worship, where background noise can interfere with speech comprehension. ALDs are designed to enhance the clarity of speech, minimize noise, and improve the listening experience in challenging acoustic conditions. While hearing aids or cochlear implants provide general amplification, ALDs serve as supplementary tools for enhancing specific listening situations.

Examples of Devices to Improve Hearing :

- **Hearing Aids:** These are small electronic devices worn in or behind the ear to amplify sound. Hearing aids are typically prescribed based on a person's specific hearing loss and are customizable to fit the individual's hearing needs.
- **Cochlear Implants:** These are surgically implanted devices that bypass damaged portions of the ear and directly stimulate the auditory nerve. They are designed for individuals with severe to profound sensorineural hearing loss and are typically used when hearing aids are no longer effective.
- **Bone-Anchored Hearing Aids (BAHA):** These devices are designed for individuals with conductive hearing loss, which occurs when sound cannot pass efficiently through the outer or middle ear. BAHA uses a small device implanted in the skull to transmit sound vibrations directly to the inner ear.

Assistive Listening Devices (ALDs):



a. FM Systems (Frequency Modulation Systems)

FM systems are one of the most widely used ALDs, especially in classroom or lecture settings. They consist of two main components:

- **Microphone (transmitter):** Worn by the speaker (e.g., teacher or lecturer).
- **Receiver:** Worn by the listener (either as a small device or integrated into a hearing aid or cochlear implant).
- FM systems use **radio frequencies** to transmit the speaker's voice directly to the listener's receiver, ensuring the sound is delivered clearly and without interference from background noise. The range and quality of the sound are significantly improved in environments such as classrooms, theaters, and public speaking venues.

b. Infrared Systems


- Infrared (IR) systems use **infrared light** to transmit sound signals from a speaker to a receiver. The speaker wears a microphone, and the sound is converted into an infrared signal, which is picked up by a receiver worn by the listener.
- **Advantages:** Infrared systems are secure and provide excellent sound quality without interference from other devices. They are often used in movie theaters, churches, and conference rooms.
- **Limitations:** The receiver and transmitter must be within the line of sight of each other for the infrared signal to be effective, which can limit the range in large venues.



Assistive Listening Devices (ALDs):

c. Induction Loop Systems

- Induction loop systems use **electromagnetic fields** to transmit sound directly to hearing aids equipped with telecoil (T-coil). These systems are installed in venues like churches, theaters, airports, and other public spaces.
- **How it works:** An induction loop (wire) is installed around the perimeter of a room or space. When a sound source (e.g., a microphone) is connected to the loop, it creates an electromagnetic signal that is picked up by hearing aids with a T-coil. The listener then hears the amplified sound directly through their hearing aids, without interference from background noise.
- **Advantages:** Induction loops are effective in public spaces and help individuals with hearing aids or cochlear implants. They provide a seamless listening experience.
- **Limitations:** The user must be within the loop's range to hear the signal, and T-coil compatibility is required in hearing aids.



d. Personal Amplifiers


- Personal amplifiers are small, portable devices that help amplify sound in one-on-one or small group settings. They consist of a microphone and a small receiver or speaker that the listener uses to hear the amplified sound.
- **Use case:** Personal amplifiers are ideal for conversations in noisy environments or small group settings where hearing aids may not be effective.
- **Advantages:** These devices are compact, user-friendly, and useful in personal interactions. They can be used without the need for a hearing aid or cochlear implant.
- **Limitations:** Personal amplifiers are generally not designed for larger spaces or more complex listening environments.



Assistive Listening Devices (ALDs):

e. Bluetooth Hearing Aid Compatibility

- Many modern hearing aids and ALDs are now **Bluetooth-enabled**. This allows users to connect their hearing aids directly to various devices, such as smartphones, TVs, computers, and other Bluetooth-compatible devices, for a clearer and more direct transmission of sound.
- **Use case:** Bluetooth functionality is particularly beneficial for people who need to stream audio from smartphones or televisions directly into their hearing aids without background noise.
- **Advantages:** Bluetooth connectivity eliminates the need for extra devices and cables, making it easier to access audio from different sources.



Applications

- Classrooms and Educational Settings:** FM systems are widely used in schools to help students with hearing loss hear the teacher's voice over background noise. They help students participate more effectively in lessons and discussions.
- Public Venues:** ALDs such as induction loop systems and infrared systems are commonly found in theaters, churches, airports, and conference rooms to improve speech clarity for people with hearing impairments.
- Personal Conversations:** Personal amplifiers help individuals with hearing loss engage in one-on-one conversations, especially in noisy environments, by amplifying the speaker's voice.
- Television and Media:** Bluetooth-enabled hearing aids and ALDs allow users to stream sound directly from their television, smartphone, or computer, improving their listening experience without external noise distractions.
- Places of Worship and Ceremonies:** Induction loop systems and FM systems are commonly used in places of worship and for ceremonies to help individuals with hearing loss participate in services, lectures, or events.

User Groups of ALDs:

ALDs are beneficial for a variety of user groups, including:

- Individuals with Hearing Loss:** People with hearing impairments who have difficulty understanding speech in noisy or large environments. ALDs are particularly useful for those who wear hearing aids or cochlear implants but still face challenges in certain situations.
- People without Hearing Loss:** ALDs can also be useful for individuals who do not have hearing loss but may need enhanced sound clarity in environments with significant background noise.
- Public Institutions and Businesses:** Schools, theaters, conference centers, and places of worship use ALDs to ensure accessibility for individuals with hearing impairments and to comply with accessibility laws.